





THE HIGH-END MICROSCOPE AT THE PRICE OF THE MIDDLE ONE




Just few minutes with SmartOPTIC will inspire you with confidence!

We are proud to introduce a SmartOPTIC medical microscope dedicated for demanding users.

Solutions applied in this microscope provide easy mobility and great comfort of work.

With superior optical quality, SmartOPTIC offers high optical resolution, accuracy and focal range - all that makes SmartOPTIC an indispensable tool in modern medical practice.

SmartOPTIC - choice of demanding professionals



WHEN HIGH EXPERTISE AND ACCURACY COUNT...



Weight Balance System

Well balanced arm provides smooth movement and precise positioning. Wide adjustment range allows for good balance also after auxiliary equipment installation.

LED Light Source

If you respect economy and very bright light with a colour temperature of 6000 Kelvins, near sunlight colour, then this Light Source is for you.

The 80 Watt diode is rated at 50000 hours.

Xenon Light Source

Xenon Light Source will suit the needs of professionals demanding the highest intensity of light.

It provides a very bright (105 kLux) and white light with a colour temperature of 6000 Kelvin, near sunlight in colour. It enables work in a natural light environment. The 180 Watt bulb is rated at 500 hours.

Its intensity is controlled by a shutter that preserves colour temperature at all brightness levels.


Halogen Light Source

Halogen Light Source has two independent channels 150 Watt each. Light with a colour temperature of 5400 Kelvins is an irreplaceable illumination for most of clinical cases.

Reliable support system

Heavy but mobile base with rubberized wheels for firm yet easy movement with breaks for fixing position.


See other mounting options on page 10.



be adjusted easily by pressing the axially located button.

200, 300 or 400 mm.

THE MORE YOU CAN SEE, THE MORE YOU CAN DO



Dedicated PROCAM digital camera



Ultra compact design with dedicated beam spliter. Auto white balance. View frame delimitation. Image capture button.

Vision track with digital camera


Basic video recording and high - resolution still images for documentation, archive and educational purposes.

Customized adaptors for SLRs and most of compact cameras available on request.



Digital camera recommended for SmartOPTIC microscopes.

Light-weighted, powerful camera with HD recording, 14MP sensor, Live View output and resonable price.



High resolution preview and recording with your personal video camera recorder e.g. SONY Handycam.













Capturing system for images, video and sound, specially designed for use in the dentist's practice.

Small, portable box that can be used with a laptop or desktop computer.



Combination of vision track equipment

Two socketed beam spliter can be armed with any combination of the above devices that can work together at the same time.

SMART ACCESSORIES VARIETY OF ENHANCEMENTS



Binocular Extender

Applied together with an inclinable binocular allows to maintain a comfortable upright position. reducing muscle strain.

Crucial ergonomy enhancement.



Binocular Rotation Ring

Provides additional comfort. Even when the head of the microscope is angled the ring gives you the possibility to maintain your upright position.



Lens cap Protecting the objective lens especially while grinding.





Filters


Green- filter to work with patients'

Orange - filter prevents the premature curing of the composite materials.

Laser filter for safe work with laser devices.

Polarizing filter - easy identification of tissue changes.





Autoclavable handle covers

Preventing the infections, making your work safer for you and your patients. Covers meet high hygienic standards.





Lenses

You can chose among variety of lenses depending on your work specification: f-200, f-250, f-300, f-400



Ceiling mount Ideal positioning for quick

Ideal positioning for quick access. Mount for ceiling height above 200 cm.







SmarloPIIC S

Table mount

For trainings but also can be useful if the stable and solid work surface is located near the operation unit.

Mobile floor stand

Construction on wheels allows the movement between sites.



MAGNIFICATION TABLE

Objective Lens		200				250				300				400			
eyepiece	position of button	magnification		field		magnification		field		magnification		field		magnification		field	
		f135	f170	f135	f170	f135	f170	f135	f170	f135	f170	f135	f170	f135	f170	f135	f170
10x/16mm	1 (0.4x)	2,70	3,40	59,26	47,06	2,16	2,72	74,07	58,82	1,80	2,27	88,89	70,48	1,35	1,70	118,52	94,12
	2 (0.6x)	4,05	5,10	39,51	31,37	3,24	4,08	49,38	39,22	2,70	3,40	59,26	47,06	2,02	2,55	79,21	62,75
	3 (1x)	6,75	8,50	23,70	18,82	5,40	6,80	29,63	23,53	4,50	5,66	35,56	28,27	3,37	4,25	47,48	37,65
	4 (1.6x)	10,80	13,60	14,81	11,76	8,64	10,88	18,52	14,71	7,20	9,06	22,22	17,66	5,40	6,80	29,63	23,53
	5 (2.5x)	16,88	21,25	9,48	7,53	13,50	17,00	11,85	9,41	11,25	14,17	14,22	11,29	8,42	10,63	19,01	15,06
	1 (0.4x)	3.37	4,25	59,35	47.06	2.70	3,40	74,07	58.82	2.25	2.83	88.89	70.67	1,69	2,13	118,34	93.90
12.5x/16mm	2 (0.6x)	5,06	6,38	39,53	31,35	4,05	5,10	49,38	39,22	3,38	4,25	59,17	47,06	2,53	3,20	79,05	62,50
	3 (1x)	8,44	10,63	23,70	18,81	6,75	8,50	29,63	23,53	5,63	7,08	35,52	28,25	4,20	5,30	47,62	37,74
	4 (1.6x)	13,50	17,00	14,81	11,76	10,80	13,60	18,52	14,71	9,00	11,32	22,22	17,67	6,74	8,50	29,67	23,53
	5 (2.5x)	21,08	26,58	9,49	7.52	16,88	21,25	11,85	9,41	14.08	17,71	14,20	11.29	10,54	13,33	18,97	15,00



SYSTEM COMPONENTS

Binoculars:

- fixed angle 45° fixed angle 90° adjustable angle 0-195°
- all with diopter adjustment

Objective lenses:

• working distance 200 mm • 250 mm • 300mm • 400mm

Illumination:

- halogen light source 3000K 37klux LED light source 6000K 50klux
- xenon light source 6000K 105 klux high performance fiber optic cable

Filters:

■ green ■ orange ■ laser filter ■ polarizing filter

Mounting options:

• mobile floor stand • wall mount • ceiling mount • table mount

Microscope arm:

standard 60 cm extended 90 cm 28 cm additional arm extension

Vision track:

- dedicated digital camera with beamsplitter
 beamsplitter
- 2-channel beamsplitter Carr binocular extender with 2-channel beamsplitter
- DSLR camera adaptor
 camcorder adaptor
 Sony NEX camera adaptor
- CCD camera adaptor monitor stand

System enhancements:

- Carr binocular extender
 binocular rotation ring
 monitor stand
- printer shelf
 autoclavable handle covers
 lens cap

BASE CONFIGURATION

Binocular: with fixed angle 45° (wide angle, diopter

adjustment)

Lens: working distance 250 mm

Light-source: halogen or LED; high performance

fiber optic cable

Filter: green or orange

Varsion: mobile floor stand

Arm length: standard 60 cm

We are always ready to **meet our clients expectations** by learning their needs and presenting the opportunities that SmartOPTIC has to offer. Our consultants will help you to **choose the best option** so you could perform your work as you have always desired.







TRAININGS OF WORKING WITH A MICROSCOPE

www.smart-optic.com



DISTRIBUIT IN ROMANIA PRIN TEHNODENT POKA

TIMISOARA - ROMANIA Str. Crișul, Nr. 2, Tel: +40.256.246.114 Mobil: +40.755.082.077 www.poka.ro | comenzi@poka-dent.ro





